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test-it IgG FOAL

Diagnostic kit to determine the titer of immunoglobulines type G (IgG) in newborn foal.

Introduction

Foals are in hypo-γ-globulinemia at birth. It is imperative that they drink and assimilate the mare's colostrum for protection. This protection against pathogens consists largely in the foal's intestinal absorption of the maternal Ig during the first hours after birth. This is the principle of transfer of passive immunity, essential to the survival of the foal.

Test Objective

The objective of this test is to determine the titer of maternal antibodies in the newborn foal. The IgG level is measured by immunochromatography solid phase (nitrocellulose) and is based on a semi-quantitative detection in foal's whole blood, serum or plasma.

Precautions and rules for optimal use (to minimize the risk of „technical“ false positive or false negative results)

- use for *in-vitro* diagnostic purposes only,
- do not use the test if its pouch is damaged or after expiration date.
- use disposable gloves and handle the samples with care, they could potentially transmit zoonotic pathogens.
- do not ever touch the result's window.
- use within 10 minutes after opening of the pouch (test is sensitive to humidity).
- use a new loop and a new diluent for each sample.
- use only the diluent and the buffer included in the box.
- use the number of drops mentioned on the box without overloading the sample well "S".
- maintain the test device horizontally throughout the duration of the test.
- do not use hemolyzed or defiled samples.
- discard wastes in accordance with local requirements.
- discard disposable gloves and wash hands after handling.
- keep out of the reach and sight of children.

Storage


- **diagnostic kits:** between +2 °C and +30 °C. Storage at room temperature (+18 °C to +30 °C) is recommended. If the tests have been refrigerated between +2 °C and +8 °C expose them to room temperature 15 minutes before use (unopened foil).
- **samples:** At room temperature, they should be tested within 4 hours after collection. If testing is delayed, samples should be kept refrigerated between +2 °C and +8 °C for up to 4 days. In this case just expose them to room temperature 15 minutes before testing. For prolonged storage, it is recommended to freeze the sample (serum or plasma only) at -20 °C. In this case, just defrost it before use.

Further examination

The interpretation of all diagnostic tests should be based on history, clinical and epidemiological context of the tested animal. These have priority regardless of the test result.

Description of the test device

Lateral Flow Immunoassay (LFI): mechanism

Principles of the test

- Test-it IgG FOAL is designed to detect the titer of IgG in newborn foal's whole blood, serum or plasma. After conducting the test:
- The "foal IgG" (if present in the sample) migrate along the nitrocellulose membrane and bind to "IgG goat-anti IgG foal" + colloidal gold particles, located in the conjugate pad (CP) → complex formation. In turn, this complex migrates on the nitrocellulose and binds to other antibodies present in the test line (T) → formation of a "direct sandwich" type complex. If the "foal IgG" are not present in the sample, nothing happens on the line T.
 - On the control line (C), a complex is formed between "Ig control + colloidal gold particles" present in CP and "control Ig" present on this C line. This immunochromatographic reaction must occur, regardless of whether the sample contains "foal IgG" or not. The coloration of this C line is a witness to the interpretation of the coloration of the T line (see interpretation of results).

Characteristics and benefits of the test-it

- display box
- rapid diagnosis: immunochromatography in one-step (LFI)
- result exclusively between 5 and 10 minutes
- storage and use at room temperature
- 5 tests per box without expensive equipment required
- high quality material, entirely produced and controlled at the same manufacturing site, which increases the sensitivity and specificity of the test.

Contents of the display box

- 5 test devices
- 5 disposable loops for sample taking
- 5 diluents
- 1 buffer

Test procedure

Sample = whole blood, serum or plasma of the newborn foal

- 1) Carefully immerse the loop (only the tip) into the sample. Prior to the next step, visually confirm that the loop is filled.
- 2) Transfer the filled loop by gently shaking its tip in the diluent.
- 3) Mix sample and diluent by inverting this one.
- 4) Hold the diluent vertically and discard the 5 first drops from the bottle.
- 5) Dispense vertically **1 drop of the diluent into the specimen well (S)** and wait for complete absorption.
- 6) Then dispense vertically **1 drop (40 µl) of the buffer into the specimen well (S)**.

Reading of the result

- exclusively between 5 and 10 minutes.
- discard the used device after 10 minutes.

Interpretation of the result

Control line (C): this purple line should always appear in the result's window after performing the test. If this line does not appear, the test should be considered invalid. In this case, repeat test again.
 Test line (T): The titer of IgG is determined in semi-quantitative way, by comparing the intensity of T line with C line (see here under).

C line versus T line		IgG level
	C >> T	< 400 mg/dL
	C > T	400 – 800 mg/dL
	C = T	= 800 mg/dL
	C < T	> 800 mg/dL

Advised protocol: 1 or 2 tests for the monitoring of the foal

- 1st test advised 12h after birth:**
- If the indicated amount of IgG is > or = 800 mg/dl, the foal is protected.
 - If the indicated amount of IgG is < 400 mg/dl, a blood or plasma transfusion can already be considered for the foal. Intestinal absorption of colostrum is still possible until 24 hours after birth, by natural way or by gastric probe.
 - If the amount of IgG is between 400 and 800 mg/dl, measures must be taken according to a risk/benefit analysis (health aspects, clinical and environmental aspects).
- 2nd test advised 24 hours after birth:** to be performed only when the first test indicates that the amount of IgG is < 400 mg/dl or between 400 and 800 mg/dl.
- If the indicated amount of IgG is > or = 800 mg/dl, the foal is protected.
 - If the indicated amount of IgG is < 400 mg/dl, a blood or plasma transfusion must be considered for the foal.
 - If the amount of IgG is between 400 and 800 mg/dl, measures must be taken according to a risk/benefit analysis (health aspects, clinical and environmental aspects).

Liability

Possible risks due to the performance of this product are assumed by the purchaser. The manufacturer is not decision maker and is not responsible for indirect, special or consequential damages of any kind resulting from the use of this product.

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